

Flexible fiberoptic
Naso-Pharyngo-Laryngoscopes.

Advanced fiber endoscopy for ENT.



The PENTAX Medical FNL product range offers a diverse range of conventional and portable fiberoptic endoscopes for high quality and cost effective endoscopic ENT examinations.

Optical instruments of the highest quality to support diagnosis and therapy.



Maximum Flexibility

- Flexible fiberoptic endoscopes can access places that cannot be reached with a rigid endoscope.
- PENTAX Medical fiberoptic endoscope design allows light source and video camera versatility to provide compatibility with a variety of different medical devices.
- Maximum mobility can be achieved by an LED battery powered cordless solution.

Broad Field of Application

- The FNL product range offers a paediatric model, routine endoscopy models and a therapeutic model with a 2.2 mm working channel.
- Both diagnostic and therapeutic procedures are within reach with the PENTAX Medical fiberoptic endoscopes.
- Optical instruments of the highest quality to support diagnosis and therapy.

High-Quality Fiberoptic Endoscopes

Type	Naso-Pharyngo-Laryngoscopes			
	Slim	Routine		Therapeutic
Product	FNL-7RP3	FNL-10RP3	FNL-10RBS	FNL-15RP3
Insertion tube [ø mm]	2.4	3.5	3.5	4.9
Distal end [ø mm]	2.2	3.2	3.2	4.8
Instrument channel [ø mm]	-/-	-/-	-/-	2.2
Working length [mm]	300	300	300	300
Field of view [°]	75	75	85	75
Tip deflection up/down [°]	130/130	130/130	130/130	130/130
Remarks	-	-	portable	-

PENTAX Europe GmbH

Julius-Vosseler-Straße 104 . 22527 Hamburg . Germany

Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13

E-mail: info.emea@pentaxmedical.com . www.pentaxmedical.com